

# Güray Bite Raiser™

Temporary Bite Opening for Deep or Cross-Bite Cases That Call for Unobstructed Tooth Movement



**Ask The Expert**  
enisg@orthotechnology.com

## Easy Chairside Application

Simply shape and mold the Güray Bite Raiser for a custom fit, insert and ligate the proximal wings into the headgear tube and ligate the "T" Spur onto the lingual cleat.

## Secure, Yet Removable

Güray Bite Raisers, unlike Bite Turbos and Bite Builders, can be easily disconnected at the lingual cleat to check the occlusion. Simply resecure if additional treatment time is needed with no debonding concerns, saving the clinician both time and money.

## FEATURES & BENEFITS


- Inexpensive and practical
- Easy and unobstructed tooth movement
- Ensure optimal oral hygiene
- Comfortable for patient
- Easy insertion and removal - insertion takes a couple of minutes, removal takes even less

## ELIMINATED

- Patient compliance requirement
- Removable bite plates
- Laboratory work
- Restorative materials

## RESULTING IN

- Less chairtime
- Less treatment time and expense



		Bite Raiser Thickness	
		0.8mm Child	1.0mm Adult
Bite Raiser Length	4.0mm Arm	Item # 10030 (for occlusal headgear tubes)	Item # 10032 (for occlusal headgear tubes)
	6.0mm Arm	Item # 10031 (for gingival headgear tubes)	Item # 10033 (for gingival headgear tubes)

### Güray Bite Raiser (Child)

0.8mm x 4.0mm .....10030  
0.8mm x 6.0mm .....10031  
10 per pack

### Güray Bite Raiser (Adult)

1.0mm x 4.0mm.....10032  
1.0mm x 6.0mm.....10033  
10 per pack



© 2013 Ortho Technology, Inc. Güray Bite Raiser is a trademark of Ortho Technology, Inc. U.S. Patent Number 6,726,473.